



Chair Bed with Riser Recliner action

The Chair Bed is a unique adaptation of the Two Motor Riser Recliner Chair and allows for independent adjustment of the leg rest and back rest. It has been designed for users who are unable to sleep in a bed. The Rotoflex is the ultimate solution to this difficulty problem for many, but the Chair bed is a possible solution too. In normal circumstances, we would not recommend anyone sleeping in a chair but for those who can't sleep in a bed and are faced with sleeping in a chair, the Chair Bed is definitely a viable option.

Why sleep in the Chair Bed rather than a Riser Recliner?

It will recline to a flat position just as a bed. This is unique to the Theraposture Chair Bed.



The arms drop down to contribute to the sleeping surface and there are no wings on the back meaning that the user can move freely on the sleeping surface just as they would on a bed. It is possible to specify a chair bed that is up to 38" wide – wider than a single bed.

Each Chair Bed is supplied with a removable topper that is filled with memory foam. This offers pressure relieving properties and adds to comfort. The topper effectively creates an uninterrupted flat sleeping surface that covers the gaps between the back and seat and leg rest and seat. It can be removed should it need to be cleaned and is often supplied in a waterproof stretchable fabric.

All Chair Beds are made to measure following an assessment. As a result of this we offer a week long suitability guarantee ensuring the products suitability for the user.

All Chair Beds have a riser function that assists with standing transfers.

This product is a solution where getting in and out of bed is no longer possible or becoming very difficult and sleeping in a chair is the only available option. (Do look at the Rotoflex page to see if this product can assist)

We provide a free of charge assessment of your needs, wishes and risks in your own home, with absolutely no obligation to purchase. To give you peace of mind, we are members of the BHTA and conform to their Office of Fair Trading approved Code of Practice, which includes standards of behaviour and ethical guidelines for member staff operating in client's homes.

All chairs have a 5 year warranty on all actuators and electronics other than the hand control. Hand controls upholstery and mechanisms have a 2 year warranty.

Hand Control

The Chair Bed with Riser Recliner action has 4 buttons allowing you to move the back rest forwards or backwards and the leg rest up or down. Once the leg rest is all of the way down, the chair will then rise to help you achieve a standing transfer. Our hand controls are the best you can get. We believe that cutting corners to save a few pounds can result in reliability issues.

Sizes Available

Action of Chair	Capacity	Width	Height	Depth	Back Height
Chair-Bed (won't go fully flat)	20 stones	26" to 30" (overall)	16.5" to 20"	18" to 22"	27" to 32"
Chair-Bed	25 stones	28" to 38" (overall)	19" to 21"	18" to 22"	27" to 32"
Chair-Bed	35 stones	32" to 38" (overall)	19" to 21"	18" to 22"	27" to 32"



Chair Bed with Riser Recliner action

Extensive range of options:

All Theraposture Riser Recliner Chairs have a huge range of styles and options.

Back Style



Button Back with Wings



Waterfall Back with Wings



Lateral Support Back with Wings



Top Split Cushion without Wings



Split Lumber Back without Wings

Arm Style



Full Scroll

This is an option for all chairs creating a more traditional look.



Half Scroll

This is standard on all chairs unless otherwise stated, creating a contemporary look.



Light Knuckle Wood

This is an option on all chairs giving a solid easily gripped end to the arm.



Dark Knuckle Wood

This is an option on all chairs giving a solid easily gripped end to the arm.



Removable Arm

An alternative to the Drop Down arm, with the solid arm lifting away on one or both sides. A real help for wheelchair/commode side transfers.



Drop Down Arms

An easy to use single pull lever is situated below the arm, allowing the arm to simply drop down to a height level with the seat cushion. You can have a drop down arm on one or both sides of your Chair or Chair Bed.

This is particularly useful when maybe transferring from a wheelchair or to allow carer access should a hoist be required. When used on a Chair Bed, the lowered arm contributes to the overall sleeping surface.

Pressure and Comfort



Personalised seat firmness is achieved by using different foam compositions depending on the personal preference and the weight of the user.



Thergel

This static foam cushion is suitable for high risk users. It is constructed from Combustion Modified High Resilient foam with a profiled (castellated) upper surface to increase pressure distribution encasing a silicon gel sac



Coolgel

is a Proprietary Silicon Polymer designed to have the same physical characteristics as human fat tissue. The softness stays unchanged at any temperature. The polymer is solid and will not flow or leak if punctured. It is also hypoallergenic - it does not contain latex or isocyanates. It is odourless and the cellular design allows air to circulate between the cells to minimise heat and moisture build up. Coolgel is classified as suitable for high risk requirements.



Memory Foam (Visco)

can be used to manage a low risk pressure redistribution requirements as well as providing wonderful comfort by moulding to your body contour. The Visco foam is used in conjunction with a high resilience foam base.

The Double Visco heat sensitive foam provides a higher degree of pressure redistribution ideal for medium risk requirements.

Both options are built into the seat and are covered with a stretchable fabric making them discreet and not affecting the finished dimensions of the chair.

Fabrics

