

## Theracare Mobile Recliner Chair



The Theracare Mobile is a fully functioning Recliner Chair using either a Dual Motor (action c) or a Single Motor Tilt in Space (action a) mechanism. If specified with a dual motor T.I.S mechanism the Theracare Mobile can also be supplied with a riser function, the footrest used in this option would be a fold down instead of the standard slide out. It is supplied with large locking castors and with a wooden push handle and slide out footrest as standard. The Theracare Mobile, as its name suggests, can be easily moved between rooms with the occupant in the chair, making this a very versatile and comfortable home mobility chair.

If specified with the built in Battery Backup, there is no need to plug the chair in when moving between rooms.



We provide a free of charge assessment of your needs, wishes and risks in your own home, with absolutely no obligation to purchase. To give you peace of mind, we are members of the BHTA and conform to their Office of Fair Trading approved Code of Practice, which includes standards of behaviour and ethical guidelines for member staff operating in client's homes.

All chairs have a 5 year warranty on all actuators and electronics other than the hand control. Hand controls upholstery and mechanisms have a 2 year warranty.



### Hand Control

The Dual Motor Chair Riser Recliner has 4 buttons allowing you to move the back rest forwards or backwards and the leg rest up or down. Once the leg rest is all of the way down, the chair will then rise to help you achieve a standing transfer. Our hand controls are the best you can get. We believe that cutting corners to save a few pounds can result in reliability issues.

Sizes Available					
Action of Chair	Capacity	Width	Height	Depth	Back Height
<b>Theracare Mobile TIS</b> (recline only)	20 stones	18" to 20"	16" to 22"*	18" to 22"	27" to 30"
<b>Theracare Mobile 2 motor</b> (recline only)	20 stones	18" to 20"	16" to 22"*	18" to 22"	27" to 30"
<b>Theracare Mobile 2 motor TIS</b>	20 stones	18" to 20"	16" to 22"*	18" to 22"	27" to 30"

\*Seat Heights are measured from the footboard

# Mobile Recliner Chair



## Theracare Mobile Recliner Chair

### Extensive range of options:

All Theraposture Riser Recliner Chairs have a huge range of styles and options.

### Back Style



**Button Back with Wings**



**Waterfall Back with Wings**



**Lateral Support Back with Wings**



**Top Split Cushion without Wings**



**Split Lumber Back without Wings**

### Arm Style



**Full Scroll**  
This is an option for all chairs creating a more traditional look.



**Half Scroll**  
This is standard on all chairs unless otherwise stated, creating a contemporary look.



**Light Knuckle Wood**  
This is an option on all chairs giving a solid easily gripped end to the arm.



**Dark Knuckle Wood**  
This is an option on all chairs giving a solid easily gripped end to the arm.



**Removable Arm**  
An alternative to the Drop Down arm, with the solid arm lifting away on one or both sides. A real help for wheelchair/commode side transfers.



**Drop Down Arms**

An easy to use single pull lever is situated below the arm, allowing the arm to simply drop down to a height level with the seat cushion. You can have a drop down arm on one or both sides of your Chair or Chair Bed.

This is particularly useful when maybe transferring from a wheelchair or to allow carer access should a hoist be required. When used on a Chair Bed, the lowered arm contributes to the overall sleeping surface.

### Pressure and Comfort



**Personalised seat firmness**  
is achieved by using different foam compositions depending on the personal preference and the weight of the user.



**Thergel**  
This static foam cushion is suitable for high risk users. It is constructed from Combustion Modified High Resilient foam with a profiled (castellated) upper surface to increase pressure distribution encasing a silicon gel sac



**Coolgel**  
is a Proprietary Silicon Polymer designed to have the same physical characteristics as human fat tissue. The softness stays unchanged at any temperature. The polymer is solid and will not flow or leak if punctured. It is also hypoallergenic - it does not contain latex or isocyanates. It is odourless and the cellular design allows air to circulate between the cells to minimise heat and moisture build up. Coolgel is classified as suitable for high risk requirements.



**Memory Foam (Visco)**  
can be used to manage a low risk pressure redistribution requirements as well as providing wonderful comfort by moulding to your body contour. The Visco foam is used in conjunction with a high resilience foam base.

The Double Visco heat sensitive foam provides a higher degree of pressure redistribution ideal for medium risk requirements.

Both options are built into the seat and are covered with a stretchable fabric making them discreet and not affecting the finished dimensions of the chair.

### Fabrics

