



Examples of the mutiple settings and the very low access



Olaf 135 natural wood showing profiling (bottom right)  
Olaf 98 natural wood 140 x 70cm (top right)  
Olaf 135 natural wood with doors open (left)



Olaf 135 natural wood raised lying surface (left)  
Olaf 98 natural wood with doors open and entry aid (top right)  
Olaf 98 natural wood with lying surface in highest position (bottom right)

## Olaf Beds - the flexible all-rounder

... an extremely low access level and a perfect care-level brought together in one bed

With two low access options and 70cm of variable height, the Olaf answers many needs

The Olaf range is highly suitable for restless children with excessive physical mobility.

It is ideal for children who can independently pull to kneel or standing and require therapy treatments in bed.

It is suitable for children, who despite illness and disability have the physical mobility to get in and out of bed independently.

The colour design of the bed makes it more attractive and stimulates the child's visual awareness. It has a light and airy feel, with plexiglass panels so that the child can see what is going on around it.

At night a secure and cosy sleeping place, during the day this Kayser Bed can be part of the child's day space and play area.

The Olaf Bed has four doors that open the long side, but 8 doors are also an option.

The door height of 98 cm or 135 cm help prevent the child from climbing out of the bed unattended.

The 70cm lift motor allows for this low access cot to raise to a safe working height - from 35 cm to 105 cm (optional 17 cm/105 cm) upper edge of mattress - This cot is ideal where the child can independently transfer because it is so low, yet it still rises high enough to ensure carers risks are addressed.

The Kayser Bed Olaf can be set at seat or wheelchair height to enable a simple and safe transfer in and out of bed, in addition to being close to the ground.

There is sufficient clearance for hoists to be used. (The option 17 cm extreme low access is however not hoist accessible).

With a weight capacity of 150 kg a parent or therapist can also join the child in the bed.


It is possible to pad an Olaf, but if padding is required, the Hannah may be a better option as the Olaf has to have removable strips of padding at either end to cope with the variable height motors.

### Reasons to buy an Olaf

- 1 Three bed sizes: 70/140 cm, 90/170 cm and 90/200 cm
- 2 Extreme height adjustable up to 70 cm, i.e. very low access level and perfect care-level
- 3 Special door lock prevents the doors from being opened from the inside
- 4 Special door fitting prevents the doors being opened from inside
- 5 Customise with fixed, manual or electric mechanism to adjust lying surface
- 6 Sweat and saliva resistant wood treatment in accordance with DIN 71-3 Childrens Toys Safety Act
- 7 Can be used as play area during the day
- 8 Plexiglass panels allow visual access both in and out
- 9 Lying surface infinitely adjustable
- 10 Custom-made units to meet the specific requirements of the patient are possible



## Sleeping Platforms (choose your requirement)



### Standard accessories

Castors plastic 13,5 cm

### Optional accessories

8 Door access  
17 cm access  
W14 Mattress  
W17 Mattress  
W24 Mattress  
Therapy mat  
Solid wood base  
Detachable posts, (2 maximum)  
Castors plastic 16.5 cm  
Castors beech 13.5 cm  
Castors plastic 10 cm  
Roller cover safety guard  
Climb out protection - wooden rods  
Climb out protection - nylon net  
Entry Aid  
Padding  
Clear Soft Windows  
Doors with plexiglass  
Pull-up rail  
Play bead rod  
Canopy  
Pedestal

### Technical Data

Door height	Olaf 98	98 cm
	Olaf 135	135 cm
Inside/Outside dimensions		140 x 70 / 162 x 82 cm
		170 x 90 / 192 x 102 cm
		200 x 90 / 222 x 102 cm
Overall height incl. castors (13,5 cm)	Olaf 98	130 cm
	Olaf 135	167 cm
Lying height at upper edge of mattress incl. castors (13,5 cm)	Olaf 98	min. 35 cm
	Olaf 135	max. 105 cm