

Double combination bed in solid natural oak

Continue to share a bed with your partner and enjoy independent function on either side of the bed for individual comfort and needs.

Your bed is specified to your exact requirements meaning that you can have any combination of functions and sizes together in a bed that looks like part of your home.

Why not choose a beautiful solid oak bed enclosure to surround the metal mechanism that powers a Theraposture bed? There is a choice of wood finishes and colours. Also you can choose between a slatted or panelled bed surround with high or low foot ends.

We also provide an extensive range of matching solid oak bedroom furniture such as bedside cabinets, chests of drawers and wardrobes. For further examples from this impressive collection, see our accessories.



Profiling only

The profiling only action will allow you to independently raise and lower both the back rest and the leg section by simply pushing a button on the hand controller. This mechanism is often the choice of partners in a double combination bed as it allows both individuals to sit up independently together. The metal construction makes this mechanism more durable for long term use and the only choice if grab handles or side rails are needed.



Variable height and profiling

In addition to the profiling only mechanism the variable height function is of significant benefit to less able people in a number of ways. If a carer or family member assists a person in bed, the variable height function will smoothly and safely raise the occupant to a safe working height, thereby reducing the risk of a back care injury.

Additionally, this function will help a person to achieve a standing transfer by vertically raising them to standing from sitting in a controlled and safe movement.

Freephone: 0800 834654

Email: info@theraposture.co.uk www.theraposture.co.uk



Rotoflex turning system

In addition to the variable height and profiling action, the remarkable Rotoflex turning system is a unique solution to an everyday problem – getting legs into and out of bed. This is achieved by a unique rotational action that will sequentially transfer you from a sitting to a sleeping position with the touch of just one button.


In many cases, the Rotoflex truly delivers independence as it is a real alternative to hoisting and carer assistance. In this situation the bed can pay for itself many times over as the cost of care in your home or in a residential home is so high.

The Rotoflex can be made to rotate to either side of the bed to ensure that it is compatible with your room situation.



Technical options – adjustable bed in solid natural oak

Electronic profiling	Yes	Double combination	Yes
Electronic variable height	Yes	Separable halves	No
Rotoflex	Yes	Single bed	Yes
2'6" wide	Yes	Compatible with cot sides	Yes
3' wide	Yes	Suitable for castors	No
3'6" wide	Yes	Suitable for static feet	Yes
4' wide	Yes	Standard length 6'6"	Yes
4'6" wide	Yes	Shorter length available	Yes
5' wide	Yes	Longer length available	Yes



Finishes

1. Natural Oiled
2. Natural Lacquered
3. Smoked
4. Milan
5. Lime
6. Olive

E&OE

Freephone: 0800 834654

Email: info@theraposture.co.uk www.theraposture.co.uk

Theraposture Limited, Kingdom Avenue, Northacre Industrial Park, Westbury, Wiltshire. BA13 4WE
Open: Monday – Friday, 8:30am–5:30pm. Answerphone messages can be left at all other times

